One-Fil® Line

Calcium Silicate based Bioceramic Materials



☑ Beyond Simple Sealing, Toward True Healing of the Tissue

**Simple! Fast!
Predictable!
Stress-Free!**

MEDICLUS CO., LTD.

MEDICLUS aims to provide a significantly faster, healthier, and more convenient dental care system

Seoul Office | +82-10-8535-2877
No. 206, 464, Cheongpa-ro, Jung-gu, Seoul 04510, Korea
Fax : +82-2-6499-1006 E-mail : sales@mdclus.com

Head Office | Research & development center
No. 1210, 134, Gongdan-ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do 28576, Korea

MEDICLUS CO., LTD.

www.mediclus.co.kr

MEDICLUS ENDODONTIC SOLUTION

ISOLATION
Dentto-Dam™

EDTA Irrigation
Endo@ Prep™ Gel

OBTURATION
One-Fil®
One-Fil® Putty
One-Fil® Putty Injectable

ETCHING
Any-Etch®
Any-Etch® Classic

BONDING
Hi-Bond® Universal

BULK-RESIN
Any-Com® Bulk

CORE BUILD-UP
Any-Core®

Line

MEDICLUS CO., LTD.

Hi-Bond® Universal

Multifunctional Universal Adhesive



☑ Enhanced Adhesion and Desensitizing Effect with 10-MDP and MBG



**Direct or Indirect!
Any Etching
Technique!**

Endodontic Solution

MEDICLUS CO., LTD.

1 ISOLATION

Dentto-Dam™


Light-cured liquid dam for precise isolation

Protect Completely, Remove in One Piece



Advantages

- Gingival protection during whitening treatment
- Gingival barrier against etchant, blasting powder, and other irritants
- Auxiliary sealing with rubber dam application
- Temporary wall formation during endodontic treatment



Isolation Rubber Dam Ⅰ | Isolation Rubber Dam Ⅱ | Endo-Wall



CONTENTS

402 Dentto-Dam 1.2ml x 5ea, Disposable tips

2 EDTA Irrigation

Endo@Prep™ Gel


EDTA Root Canal Cleaning and Preparation

Root Canal Irrigation, Cleaner Than Ever



Advantages

- Removal of inorganic substances through chelation ($EDTA\ Ca^{2+}$)
- Easy application with lubricating action
- Provides a clean, smear-free canal wall



<References> Jeonbuk National University School of Dentist Conservative Dentistry Prof. Kyung-San Min's Treatise



CONTENTS FDA

501-1 Endo@Prep EDTA Gel 9g x 2ea

3 OBTURATION


One-Fil® Line


Calcium Silicate based Bioceramic Materials

One-Fil®
One-Fil® Putty
One-Fil® Putty Injectable

DA DENTAL ADVISOR ★★★★½

FROM SEALING TO HEALING





Advantages

- Excellent Flowability & Wettability for Complex Root Canal Systems
- Long-Term Sealing Through Excellent Volumetric Stability
- Bioactive Healing Environment
- Reliable Curing & Long-Term Sealing Performance
- Simplified Root Canal Procedure
- Less Waste per Application



CONTENTS FDA

508-1 One-Fil 2g x 1ea



CONTENTS FDA

509-1 One-Fil Putty Injectable 0.5g x 4ea



CONTENTS FDA

509-2 One-Fil Putty 0.25g x 4ea

4 ETCHING

Any-Etch®
Any-Etch® Classic

37% Phosphoric acid

Etching Gel for Strong, Reliable Bonding



Advantages

For pretreatment of teeth and restorations in adhesive procedures : Composite resin restoration, Adhesive bonding of all-ceramic and composite resin indirect restorations, Splint, Orthodontic appliances, Tooth jewelry, Sealants

- Select the Ideal Viscosity



Any-Etch® | Any-Etch® Classic



CONTENTS FDA

304-1 Any-Etch 3ml x 5ea

305-1 Any-Etch Classic 3ml x 5ea

5 BONDING

Hi-Bond® Universal

DA DENTAL ADVISOR ★★★★½

Multifunctional Universal Adhesive

One Solution for Bonding and Desensitizing



Advantages

- Excellent Wettability for Effective Dentinal Tubule Penetration
- Excellent Compatibility with Dual-Cure Resin Materials
- Easy Handling with One-Bottle System and Drop Control Design
- Biological Safety Supported by Cell Compatibility
- MBG (Mesoporous Bioactive Glass) - Based Bioactivity
- Tooth Remineralization Effect



CONTENTS FDA PATENTED

303-2 Hi-Bond Universal 5ml x 1ea

6 BULK-RESIN

Any-Com® Bulk


Fluoride-releasing Bulk-fill Flowable Resin

Fast Filling, Reliable Protection



Advantages

- Save time
- One layer curing up to 4mm
- Simple and Easy Technique
- Excellent Fluoride Release for Pediatric Applications
- Excellent Radiopacity





CONTENTS SHADE | A1, A2, A3, TL, Blue FDA

#307-11 Any-Com Bulk A1 2g x 1ea

7 CORE BUILD-UP

Any-Core®

Dual-cured Core Build-up Composite Resin

Strong Protection, Gentle Care



Advantages

- Freedom to choose for efficient esthetics
  - Three different shades to meet clinical needs : A3, White, and Blue.
  - Two curing modes : Light-Cure and Self-Cure
- Robust core build-up
  - Convenient and fast clinical process in conjunction with Hi-Bond Universal. Proven core build-up system provided by a trusted manufacturer



CONTENTS SHADE | A3, White, Blue FDA

801-1 Any-Core Kit A3 9g x 2ea